
140119-8






- English
- Deutsch
- Français
- Italiano

QuikRead go[®]

CRP+Hb

Assay Procedure • Testverfahren

8

| Measure | | RESULT |
|--|----------|---------------------------------------|
| CRP: | | 20 mg/l |
| Hb: | | 125 g/l |
| Patient ID: | XXXXXXXX | Measurement time: 2012-06-11 12:19 |
| Test: | CRP + Hb | Result info |
| i Choose Result info to view result information. Remove cuvette to perform a new measurement. | | |
| Exit | Print | New measurement |

1 Intended use

The QuikRead go CRP+Hb test is intended for quantitative determination of CRP (C-reactive protein) in whole blood, serum and plasma and for quantitative determination of hemoglobin (Hb) in whole blood using the QuikRead go instrument. For *in vitro* diagnostic use only.

2 Summary and explanation of the test

CRP

CRP is an acute phase protein present in low concentrations in healthy individuals¹. Any pathological condition associated with invasive bacterial infection, inflammation or tissue destruction is accompanied by elevation of the CRP level in the patient's serum. The rise in CRP levels is rapid, and increased levels can be detected within 6 to 12 hours of the onset of the inflammatory process².

Quantitative measurement of the CRP concentration has been reported to be a sensitive indicator of the efficacy of antimicrobial therapy and the course of bacterial infections, as well as an effective tool in controlling and monitoring postoperative infections²⁻⁶.

Hemoglobin

Hemoglobin is the iron-containing protein in red blood cells that carries oxygen. Determining its concentration is a common procedure in both primary health care and acute care. A decrease in the hemoglobin concentration may follow from an iron deficiency or loss of blood. An elevated concentration of hemoglobin in blood may indicate a low oxygen supply, e.g. in those living in mountain regions or in smokers.

3 Principles of the procedure

CRP

The CRP measurement of the QuikRead go CRP+Hb test is turbidimetric and it is based on the agglutination reaction between microparticles coated with anti-human CRP F(ab)₂ fragments. The CRP present in the sample reacts with the microparticles, and the resultant change in the turbidity of the solution is measured by the QuikRead go instrument. QuikRead go CRP correlates well with other immunoturbidimetric methods.

The buffer haemolyses the blood cells in a cuvette, and the agglutination reaction is measured after adding the reagent. Test calibration data is in the bar code label of the cuvette and it is read by the QuikRead go instrument before the test is started. The CRP value is automatically corrected according to the sample's hematocrit level. Hematocrit results are not displayed but used in calculations.

Hemoglobin (Hb)

The hemoglobin measurement using the QuikRead go CRP+Hb test is based on measuring oxyhemoglobin photometrically at two wavelengths. The buffer haemolyses the sample's red blood cells in the cuvette, whereupon the absorption of released hemoglobin can be measured. The instrument corrects the effect of lipids or other turbidity-inducing factors on the results. The bar code on the cuvette contains the required parameters for calculating the hemoglobin concentration, and the QuikRead go instrument displays the obtained hemoglobin results.

4 Reagents

Content of the kit

| Component name and origin | Symbol | QuikRead go [®] CRP+Hb Cat. no. 140068, 50 tests |
|--------------------------------|-----------|--|
| CRP reagent caps ORIG SHP | REAG CPS | 2 x 25 |
| Cuvettes prefilled with buffer | BUF | 2 x 25 x 1 ml |
| Capillaries (20 µl) | CAPIL HEP | 50 |
| Plungers | PLUN | 50 |
| Instructions for use | | |

Reagents contain preservatives, please see Section 5 "Warnings and precautions".

Reagent storage

| Reagent | Storage at 2...8 °C | Storage at 18...25 °C |
|---|---|--|
| CRP Reagent caps (in opened and unopened aluminium tube) | Until the expiry date of the reagent caps | 24 h per day – 1 month 7.5 h per day – 3 months |
| Prefilled cuvettes in unopened foil pouches | Until the expiry date of the prefilled cuvettes | Until the expiry date of the prefilled cuvettes |
| Unopened prefilled cuvettes without the foil pouch | 6 months | 3 months |
| Opened prefilled cuvette | 2 hours | 2 hours |

Mark the date of opening the foil pouch on the cuvette rack.

Reagent preparation and storage conditions

All reagents are ready to use. Keep the CRP reagent caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

Reagent deterioration

The product should be used only if the buffer volume in the cuvette is correct. Check that the liquid surface is between the two lines marked on the cuvette. Do not use a cuvette with visible dirt in the buffer.

5 Warnings and precautions

Health and safety information

- For *in vitro* diagnostic use only.
- Do not smoke, eat or drink in premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of soap and water.
- All patient samples and controls should be handled as potentially infectious material.
- The buffer contains 0.004% reaction mass of: 5-chloro-2-methyl-2H-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) (Skin Sens 1, Aquatic Chronic 3), and < 0.1% sodium azide. May cause an allergic skin reaction (H317). Harmful to aquatic life with long lasting effects (H412). Avoid breathing vapours (P261). Avoid release to the environment (P273). Wear protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents according to national and local law (P501).
- Lyophilised reagents inside a reagent cap contain < 1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501). Contact with acids liberates very toxic gas (EUH032).
- Reconstituted and liquid reagents contain < 0.1% sodium azide, which is not considered a harmful concentration. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Disposal: see Chapter 14.

Warning

Analytical precautions

- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- The QuikRead go CRP+Hb kit is intended solely for use with the QuikRead go instrument.
- Do not mix components from different lot numbers or different tests. The components are disposable; never reuse components already used for performing a test.
- When opening a kit for the first time, ensure that the foil pouches protecting the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes inside it. In addition, always ensure before using an individual cuvette that its foil covering is intact.
- Do not touch the clear flat surfaces on the lower part of the cuvette (optical part). Discard any cuvettes with fingerprints.
- QuikRead CRP Reagent Caps are colour-coded blue to distinguish them from other QuikRead analytes.

- Keep the QuikRead CRP Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.
- Do not splash liquid into the measuring well of the instrument.

6 Sample collection and preparation

Sample material, sample volume and collection

| Sample material | Sample volume | Obtainable result | Sample collection |
|--------------------------------|------------------|--------------------|---|
| Fingertip sample (whole blood) | 20 µl | CRP and hemoglobin | Prick the clean and dry finger with a lancet and discard the first drop. Wipe the finger and collect 20 µl of blood from the second drop in a heparinised capillary. |
| Anticoagulated whole blood | 20 µl | CRP and hemoglobin | Use a venous blood sample collected in a tube containing heparin or EDTA. Mix whole blood by inverting the tube several times and take 20 µl of the sample with a capillary or pipette. |
| Plasma | 20 µl or 12 µl * | CRP | Use EDTA/heparin plasma. Avoid clearly haemolysed samples. Mild haemolysis will not affect the test result. |
| Serum | 20 µl or 12 µl * | CRP | Avoid clearly haemolysed samples. Mild haemolysis will not affect the test result. |

* If you wish to use the 12 µl sample volume and achieve a wider measurement range, change the volume setting of the QuikRead go instrument. Consult the instrument instructions for use.

Sample dilution

| Sample material | Instructions |
|-----------------|--|
| Whole blood | Do not dilute whole blood samples. |
| Plasma/serum | Plasma or serum samples may be diluted with a 0.9% NaCl solution before adding the sample to the cuvette. The recommended dilution ratio is 1+3. Multiply the obtained result by four (4). |

Sample storage

| Sample material | Short-term storage | Long-term storage |
|---|--|---|
| Fingertip blood (whole blood) | Max. 15 minutes in the heparinised capillary | Do not store |
| Anticoagulated whole blood | 2...8 °C for 3 days | Separate plasma and store at below -20 °C |
| Plasma | 2...8 °C for 7 days | Below -20 °C |
| Serum | 2...8 °C for 7 days | Below -20 °C |
| Sample (whole blood, plasma, serum) in buffer | 18...25 °C for 2 hours | Do not store |

Allow the samples to reach room temperature (18...25 °C) before testing. Frozen samples must be thawed out, allowed to reach room temperature and mixed carefully prior to testing. If the thawed plasma or serum samples contain clots, they must be centrifuged. Samples should not be repeatedly frozen and thawed.

7 Procedure

Materials required but not supplied

| Material | Cat. no. | Additional info |
|--|----------------|---------------------------------|
| QuikRead go® Instrument with software version 4.1 or newer | 133893, 149915 | |
| QuikRead go® CRP Control | 153764 | Recommended for quality control |
| QuikRead go® CRP Control High | 153763 | Recommended for quality control |
| QuikRead go® Hb Control | 141154 | Recommended for quality control |
| Fingertip lancets | | |

Assay procedure

Open the foil pouch protecting the cuvette rack and mark the opening date on the cuvette rack.

Sample collection (Fig. 1–5)

Glass capillaries (20 µl) and plungers are recommended for adding the sample into a prefilled cuvette.

- 1 A prefilled cuvette must reach room temperature (18...25 °C) before use. This will take 15 minutes for an individual refrigerated (2...8 °C) prefilled cuvette. Do not touch the clear flat surfaces in the lower part of the cuvette (the optical part). Remove the foil cover from the cuvette. Beware of splashing the liquid. Condensed buffer on the foil cover has no effect on results. The test must be performed within two (2) hours from opening the cuvette. The liquid surface should lie between the two lines marked on the cuvette.
- 2 Insert the plunger into the capillary at the end with the blue stripe.
- 3 Fill the capillary with the sample up to the white stopper (20 µl). Ensure that there are no air bubbles in the capillary. Use a soft cloth or tissue to wipe away any residual sample from the outside of the capillary.
- 4 Place the sample capillary into the buffer solution in the cuvette and dispense by pressing down the plunger. Make sure that the capillary is completely empty.
- 5 Close the cuvette tightly with a CRP reagent cap. Do not press down the blue inner part of the reagent cap. After the sample has been added to the buffer, the solution is stable for two (2) hours. Keep the cuvette in an upright position, do not shake.

Analysing the sample (Fig. 6–8)

For full instructions, see QuikRead go instrument instructions for use. The display will guide you through the test.

- 6 On the QuikRead go instrument display, select **Measure**.
- 7 Place the cuvette in the instrument with the barcode facing you, as shown in Figure 7. The instrument recognises the sample type. The display shows how the measuring progresses. The instrument first measures the sample blank and then the CRP and Hb concentrations within one minute.
- 8 The result appears on the display when the measurement is completed and the cuvette automatically rises from the measurement well.

Analysing the control (Fig. 7–8)

Precise instructions on the use of controls are provided in the package insert of the control vials. Analyse the control samples as you would clinical samples, but select **Quality Control** on the instrument display (Fig. 7–8). Select the control that matches the test you are using (CRP or Hb). The result will be saved as a control measurement.

8 Quality control

Regular use of QuikRead go CRP Control (Cat. no. 153764) and/or QuikRead go CRP Control High (Cat. no. 153763) and QuikRead go Hb Control (Cat. no. 141154) is recommended.

QuikRead CRP controls: All QuikRead CRP controls are ready-to-use. The control values have been determined for the QuikRead go instrument, and the control is applied using the same procedure as for a plasma/serum sample. The sample volume is 20 µl. If you wish to use the 12 µl sample volume, change the volume setting (plasma/serum) of the QuikRead go instrument. Consult the instrument instructions for use. The control is stable for 15 minutes in the cuvette.

Other commercially available CRP controls: Handle and use the control as instructed and measure it as you would the QuikRead CRP controls. Acceptance limits for the control value must be defined with the QuikRead go instrument. The blanking process might not succeed if you use a control that contains artificial red blood cells as these might not haemolyse normally.

QuikRead go Hb control: QuikRead go Hb Control is a ready-to-use control. The control values have been determined for the QuikRead go instrument, and the control is used as whole blood sample. The sample volume is 20 µl.

Other commercially available Hb controls: Handle and use controls as instructed. Acceptance limits for the control value must be defined with the QuikRead go instrument. The blanking process might not succeed if you use a control that contains artificial red blood cells as these might not haemolyse normally. Controls in which the hemoglobin is present in a form other than oxyhemoglobin, or which contain compounds which may modify hemoglobin, might yield results that are not aligned with those obtained by other methods.

9 Results' interpretation

CRP

Elevation of CRP value is unspecific and the results should be interpreted in the light of other clinical findings.

| CRP test result | Interpretation of result ⁷ |
|-----------------|---|
| < 10 mg/l | Excludes many acute inflammatory diseases, but does not specifically exclude inflammatory processes. |
| 10–50 mg/l | Elevated concentrations found in acute disease that occur in the presence of a slight to moderate inflammatory process. |
| > 50 mg/l | Indicates high and extensive inflammatory activity. |

Hemoglobin

Hemoglobin results can only be obtained from whole blood samples.

| Reference group | Reference values ¹⁰ |
|----------------------|---|
| Women | 120–150 g/l (12.0–15.0 g/dl, 7.4–9.3 mmol/l) |
| Men | 130–170 g/l (13.0–17.0 g/dl, 8.1–10.6 mmol/l) |
| Children 1–12 years* | 110–155 g/l (11.0–15.5 g/dl, 6.8–9.6 mmol/l) |

* Children's hemoglobin values increase progressively reaching gradually adult levels.

Result levels are affected by patients' diet and the elevation from sea level of their place of residence, so it is recommended that each laboratory determines their own reference values for hemoglobin. When diagnosing patients, the patient's overall clinical status should always be taken into account.

10 Limitations of the procedure

Assay procedures other than those specified in these instructions may yield questionable results. Some substances may interfere with the test results; please see Section 12. "Performance Characteristics".

Test results should never be used alone, without a complete clinical evaluation, when making a diagnosis. Intra-individual variations in CRP are significant and should be taken into account – for example, by means of serial measurements – when interpreting the values.

11 Expected values

The assay reference range limit was determined by using 143 apparently healthy adults (59 males and 84 females) with age 19–65 according to CLSI EP28-A3C guideline. Based on the results the 95% reference limit was concluded to be ≤ 5 mg/l for each sample type. Results are based on a non parametric approach.

However, each laboratory is recommended to establish a range of normal values for their population in their region.

12 Performance characteristics

CRP

Method comparison

Patient plasma samples were measured using two clinical laboratory methods, and the QuikRead go CRP method. Summary of correlation studies is presented in the table below.

| Passing-Bablok analysis | |
|---|---|
| Clinical laboratory method 1 | Clinical laboratory method 2 |
| $y = 1.00x + 1.0$ $r = 0.994$ $n = 116$ | $y = 0.93x + 1.4$ $r = 0.987$ $n = 113$ |

Whole blood vs. plasma comparability

In a comparison of 104 clinical samples, whole blood (y) and plasma samples (x) yielded comparable results. $y = 1.03x - 0.3$

Precision and reproducibility

A precision study was performed according to Clinical and Laboratory Standards Institute (CLSI) guideline EP5-A2.

| Within-run, between-day and total precision | | | | | | |
|---|---------------|----------------|-----------------|-------------------|--------------------|--------------|
| Sample material | Sample number | Number of days | Mean CRP (mg/l) | Within-run CV (%) | Between-day CV (%) | Total CV (%) |
| Whole blood | Sample 1 | 20 | 9 | 4.5 | 4.3 | 6.2 |
| | Sample 2 | 20 | 52 | 1.4 | 2.4 | 3.2 |
| | Sample 3 | 20 | 177 | 1.9 | 2.8 | 3.6 |
| Plasma | Sample 1 | 20 | 10 | 3.7 | 0.1 | 4.2 |
| | Sample 2 | 20 | 55 | 1.5 | 0.4 | 1.9 |
| | Sample 3 | 20 | 148 | 2.0 | 1.0 | 2.4 |
| Control | Sample 1 | 20 | 30 | 2.0 | 3.1 | 3.7 |

Interference

| Interfering substance | No interference found up to concentration |
|-------------------------------------|---|
| Bilirubin | 400 μ mol/l |
| Vitamin C | 200 μ mol/l |
| Triglycerides | 11.5 mmol/l |
| Cholesterol | 9.0 mmol/l |
| Rheumatoid factors (RF) | 525 IU/ml |
| Leucocytes | 145x10 ⁹ cells |
| Anticoagulants (Li-heparin or EDTA) | No interference |

Most heterophilic or anti-sheep antibodies in the samples do not interfere with the test, as the assay antibodies lack the Fc-part. In rare cases, interference by IgM myeloma protein has been observed.

Antigen excess

CRP concentrations less than 1,000 mg/l do not give falsely low results.

Hemoglobin

Method comparison

Patient whole blood samples were analysed using a reference method^{8, 9}, a point-of-care test and the QuikRead go CRP+Hb test. A summary of correlation studies is presented in the table below.

| Passing-Bablok analysis | |
|--|---|
| Reference method ^{8, 9} | Point-of-care test |
| $y = 1.07x - 7.8$ $r = 0.996$ $n = 59$ | $y = 1.01x - 1.1$ $r = 0.994$ $n = 130$ |

Precision and reproducibility

A precision study was performed according to Clinical and Laboratory Standards Institute (CLSI) guideline EP5-A2.

| Within-run, between-day and total precision | | | | | | |
|---|---------------|----------------|---------------|-------------------|--------------------|--------------|
| Sample material | Sample number | Number of days | Mean Hb (g/l) | Within-run CV (%) | Between-day CV (%) | Total CV (%) |
| Whole blood | Sample 1 | 20 | 96 | 1.7 | 1.2 | 2.2 |
| | Sample 2 | 20 | 134 | 1.7 | 1.3 | 2.1 |
| | Sample 3 | 20 | 164 | 1.5 | 0.9 | 2.2 |
| Control | Level 1 | 7 | 126 | 1.6 | 2.1 | 3.0 |

Interfering substances

| Interfering substance | No interference found up to concentration |
|---------------------------------------|---|
| Bilirubin | 200 μ mol/l |
| Vitamin C | 200 μ mol/l |
| Triglycerides | 11.5 mmol/l |
| Rheumatoid factors (RF) | 525 IU/ml |
| Leucocytes | 65 x 10 ⁹ /l |
| Thrombocytes | 1000 x 10 ⁹ /l |
| β -carotene | 3.7 μ mol/l |
| Ibuprofen | 500 mg/l |
| Salicylic acid / Acetylsalicylic acid | 650 mg/l |
| Tetracyclines | 200 mg/l |
| Urea | 5000 mg/l |
| Anticoagulants (Li-heparin or EDTA) | No interference |

Measuring ranges

QuikRead go software versions ≥ 4.1 give both the CRP and Hb results.

CRP

For the whole blood samples the measuring range is 5–200 mg/l CRP at the normal hematocrit level of 40%. If the hematocrit is lower or higher than 40%, the measuring range will change according to the following table. If the level of hematocrit is outside the range of 15–75 %, the instrument will not display the CRP result.

| Sample type | Hematocrit % | Sample volume | CRP measuring range mg/l | |
|---------------------|--------------|---------------|--------------------------|---|
| Whole blood sample | 15–19 | 20 μ l | 5–150 | In cases when the result is above or below the measuring range, the result will be displayed, for example "> 200 mg/l CRP". |
| | 20–28 | | 5–160 | |
| | 29–35 | | 5–180 | |
| | 36–41 | | 5–200 | |
| | 42–46 | | 5–220 | |
| | 47–50 | | 5–240 | |
| | 51–53 | | 5–260 | |
| | 54–56 | | 5–280 | |
| | 57–60 | | 5–300 | |
| | 61–62 | | 5–330 | |
| | 63–64 | | 5–340 | |
| | 65–66 | | 5–360 | |
| | 67–68 | | 5–390 | |
| | 69–70 | | 5–410 | |
| 71–72 | 5–440 | | | |
| 73–74 | 5–470 | | | |
| 75 | 5–510 | | | |
| Plasma/Serum sample | – | 20 μ l | 5–120 | |
| | | 12 μ l | 5–200 | |

Hemoglobin (Hb)

The measuring range of hemoglobin in whole blood samples is 50–245 g/l (5.0–24.5 g/dl, 3.1–15.2 mmol/l).

13 Traceability of test methods

CRP

The CRP assay included in the QuikRead go CRP+Hb test is standardised according to the ERM[®]-DA 472 reference material.

Hemoglobin

The hemoglobin assay in the QuikRead go CRP+Hb test is traceable to the ICHS (cyanmethemoglobin) standard 1995 and the CRM BCR-522 reference material.

14 Disposal

- Dispose of contents according to national and local law.
- All patient samples, sampling devices, controls, used tubes, caps and cuvettes should be handled and disposed of as potentially infectious material.
- Materials of the components of the QuikRead 101 and QuikRead go kits:
 - Paper: Instructions for use
 - Cardboard: Kit box including its inner parts
 - Plastic: Cuvettes, reagent caps, foil covering cuvette rack, cuvette rack, plungers, extraction vials and tubes, swabs, plunger and capillary tubes
 - Glass: Capillaries
 - Metal: Reagent cap tubes, cuvette lids, plunger and capillary tube caps
 - Not to be recycled: Lids of reagent cap tubes (several) and magnetic cards (PVC)
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

15 Troubleshooting

Error messages

Error messages displayed by the QuikRead go instrument are listed below. For more detailed information on error messages, consult the QuikRead go instrument instructions for use.

| Error messages | Corrective action |
|--|--|
| Measurement prohibited. Please check the reagent cap. | Check that the cuvette has a reagent cap and the inner blue part of the cap is not pressed down. |
| Measurement prohibited. Cuvette temperature too low. | Let the cuvette warm up to room temperature (18...25 °C). Test the same cuvette again. |
| Measurement prohibited. Cuvette temperature too high. | Let the cuvette cool down to room temperature (18...25 °C). Test the same cuvette again. |
| Test cancelled. Blank too high. | Test the same cuvette again. The blanking process has not been completed or the sample may contain interfering substances. In the latter case, the test cannot be completed. |
| Test cancelled. Unstable blank. | |
| Test cancelled. Error in reagent addition. | Perform a new test. There has been some problem during reagent addition. Make sure that the cap is properly closed. |
| No result. Hematocrit too low. | If the level of hematocrit in a whole blood sample is below the hematocrit measurement range of the instrument (see Section 12), no result is shown. |
| No result. Hematocrit too high. | If the level of hematocrit in whole blood samples is above the hematocrit measurement range of the instrument (see Section 12), no result is shown. |

Unexpected low and high results









Possible reasons for unexpected low and high results are listed in the table below.

| Problem | Possible cause | Corrective action |
|-----------------------------|---|--|
| Unexpected low result. | Too small sample volume. | Perform a new test. Make sure that the capillary is completely filled. Avoid air bubbles. |
| Unexpected high result. | Too large sample volume. | Perform a new test. Make sure that the capillary is properly filled. Wipe away any excess sample from the outside of the capillary. |
| | Too low buffer volume. | Perform a new test. Make sure that the volume is correct by checking that the liquid surface is between the two lines marked on the cuvette. |
| | The cuvette is dirty. | Perform a new test. Do not touch the clear flat surfaces in the lower part of the cuvette. |
| Unexpected low/high result. | Component from different kit lots or tests are used. | Perform a new test. Make sure that all the reagents are from the same kit lot. |
| | Incorrect reagent storage. | Perform a new test. Make sure that the reagents are stored according to the instructions for use. |
| | Incorrect sample volume: plasma/serum/ control sample | Check, that the used sample volume (20 or 12 μ l) matches the set plasma/serum sample volume on the instrument. |

References • Referenzen

- 1 Price CP, Calvin J, Walker SA, Trull A, Newman DJ, Gorman EG. Rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. *Clin Chem Lab Med* 1999, 37, 109–113.
- 2 van Leeuwen MA, van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. *Baillière's Clinical Rheumatology* 1994, 531–552.
- 3 Olaison L, Hogevik H, Alestig K. Fever C-reactive protein and other acute-phase reactants during treatment of infective endocarditis. *Arch Intern Med* 1997, 157, 885–892.
- 4 Peltola H, Unkila-Kallio L, Kallio MJT. Simplified treatment of acute staphylococcal osteomyelitis of childhood. *Pediatrics*, 1997, 99, 846–850.
- 5 Philip AGS, Mills PC. Use of C-reactive protein in minimizing antibiotic exposure: Experience with infants initially admitted to a well-baby nursery. *Pediatrics*, 2000, 106, 1–5.
- 6 Pepys MB. The acute phase response and C-reactive protein. In: Warrell DA, Cox TM, Firth JD, Benz EJ, editors. *Oxford Textbook of Medicine*, 4th ed. Oxford University Press, 2003. Vol 2, p. 150–156.
- 7 Thomas L. C-reaktives Protein (CRP) In: *Labor und Diagnose*, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft, 1992, 780-790.
- 8 International Council for Standardisation in Haematology: Expert Panel on Haemoglobinometry. Recommendations for reference method for haemoglobinometry in human blood (ICSH standard 1995) and specifications for international haemoglobinocyanide standard (4th edition). *J Clin Pathol* 1996, 49, 271-274.
- 9 NCCLS H15-A3, vol 20 No.28., Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard – Third Edition.
- 10 Dacie and Lewis. *Practical Haematology*, 11th edition, 2012.

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QuikRead go® is a registered trademark of Aidian Oy.

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




- 1 Price CP, Calvin J, Walker SA, Trull A, Newman DJ, Gorman EG. Rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. *Clin Chem Lab Med* 1999, 37, 109–113.
- 2 van Leeuwen MA, van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. *Baillière's Clinical Rheumatology* 1994, 531–552.
- 3 Olaison L, Hogevis H, Alestig K. Fever C-reactive protein and other acute-phase reactants during treatment of infective endocarditis. *Arch Intern Med* 1997, 157, 885–892.
- 4 Peltola H, Unkila-Kallio L, Kallio MJT. Simplified treatment of acute staphylococcal osteomyelitis of childhood. *Pediatrics*, 1997, 99, 846–850.
- 5 Philip AGS, Mills PC. Use of C-reactive protein in minimizing antibiotic exposure: Experience with infants initially admitted to a well-baby nursery. *Pediatrics*, 2000, 106, 1–5.
- 6 Pepys MB. The acute phase response and C-reactive protein. In: Warrell DA, Cox TM, Firth JD, Benz EJ, editors. *Oxford Textbook of Medicine*, 4th ed. Oxford University Press, 2003. Vol 2, p. 150–156.
- 7 Thomas L. C-reaktives Protein (CRP) In: *Labor und Diagnose*, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft, 1992, 780-790.
- 8 International Council for Standardisation in Haematology: Expert Panel on Haemoglobinometry. Recommendations for reference method for haemoglobinometry in human blood (ICSH standard 1995) and specifications for international haemoglobinocyanide standard (4th edition). *J Clin Pathol* 1996, 49, 271-274.
- 9 NCCLS H15-A3, vol 20 No.28., Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard – Third Edition.
- 10 Dacie and Lewis. *Practical Haematology*, 11th edition, 2012.

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