

AIDIAN



Easy, quick and accurate rapid tests.

bioLytical Laboratories INSTI®

Easy, quick and accurate rapid tests.

INSTI® tests

INSTI tests are a family of rapid *in vitro* qualitative tests for the detection of antibodies from whole blood, fingerstick blood, serum or plasma. INSTI® tests use a unique flow-through technology enabling immediate and reliable results. INSTI® tests are intended for use by trained personnel in various health care settings as assays capable of providing test results in 60 seconds.

- Specificity and sensitivity $\geq 99\%$
- Built-in IgG capture procedural control
- Lancet, pipette and alcohol swab included
- No timers required

Multiplex HIV-1/HIV-2/Syphilis Antibody Test

One sample, one minute, two results.



HIV-1/HIV-2 Antibody Test

The world's fastest HIV test



HCV Antibody Test

The first one-minute Hepatitis C Antibody Test



Contact info

Aidian Oy

aidian@aidian.eu / +358 10 309 3000

For orders: orders@aidian.eu

www.aidian.eu

Product	Product no	Order no
INSTI® HIV-1/HIV-2 Antibody test, 1 pc	90-1015	146149
INSTI® Multiplex HIV-1/HIV-2 Syphilis test, 1 pc	90-1028	147733
INSTI® HCV Antibody test, 1 pc	90-1062	155014
INSTI® HIV-1/HIV-2 Control, 20 tests	90-1036	149874
INSTI® SYPH Control, 20 tests	90-1032	154850
INSTI® HCV Test Controls, 20 tests	90-1067	155110
INSTI® MICROSAFE® pipettes, 100 pc	1050-100	155111

AIDIAN