



Simple on-site urine culturing. Uricult® Vet Model Chart

Uricult® Vet

Accurate and easy to use bacterial culture method for the diagnosis of urinary tract infections (UTI).

The Uricult Vet dipslide is for detection of microbes causing urinary tract infections in animals. It contains two culture media: CLED (green) is for detection of total bacterial count and MacConkey (red) is for detection of gram-negative microbes.

EASY TO USE Dip – Incubate – Interpret

ON-SITE TEST No need to transport the urine sample

SAFE TRANSPORTATION Maintains integrity of the sample

LONG SHELF-LIFE Can be stored at room temperature for nine months

Make your practice stand out

Uricult Vet offers veterinary clinics a fast and convenient on-site method to identify UTI, with accurate results available in 16–24 hours. This practise helps vets to avoid sending negative samples to reference laboratory and to prevent prescribing unnecessary antibiotics..

Further information: www.aidian.eu

For more detailed information, please consult the instructions for use.



Product list	Cat. No.
Uricult Vet 10 tests	06064
Cultura M 230 V	129997 / 05866
Cultura M 110/115 V	132754 / 06266
Cultura M 230 V (with a UK plug)	133192
Cultura M 12 V (with a car plug)	132697
Thermometer, accessory	129981 / 05724
Multirack, accessory	139248

It's so easy to use!



Dip, pour or drop



Incubate



Interpret

Uricult[®] Vet

Organisms on media



Colony densities *Escherichia coli* CFU/ml

